

ASPEN VISTA® CTO4 UPGRADE KIT

- CERVICAL THORACIC ORTHOSIS UPGRADE KIT
**ZERVIKAL-THORAKALE ORTHESE ERWEITERUNGSSETZ
 KIT DE ACTUALIZACIÓN DE ÓRTESIS CERVICAL TORÁCICA
 KIT DI AGGIORNAMENTO PER L'ORTESI TORACICA CERVICALE
 KIT DE MISE À NIVEAU DES ORTHÈSES CERVICO-THORACIQUES
 CERVICALE THORACALE ORTHÈSE UPGRADE KIT
 SADA PRO UPGRADE KRČNI HRUDNI ORTÉZY
 OPGRADERINGSÆT TIL HALSHVIRVEL BRYSTKASSE ORTOSE
 ЗЕР КИ АΝΑΒΑΘΜΙΣΗΣ ΑΥΧΗΝΙΚΗΣ ΘΩΡΑΚΙΚΗΣ ΟΡΘΩΣΗΣ
 CERVIKALT OPGRADERINGSSET FOR BRYSTORTOSE
 ORTEZA PIERSIOWO-SZYJNA ZESTAW ROZSZERZONY
 UPPGRADERINGSÅTS FÖR CERVIKAL THORAKAL ORTOS
 КОМПЛЕКТ ДЛЯ РАСШИРЕНИЯ ШЕЙНОГО КОРСЕТА
 KIT DE ATUALIZAÇÃO DE ÓRTESE CERVICAL TORÁCICA
 颈椎畸形器升级套件**

ASPEN VISTA MULTIPOST
BACK PANEL

ASPEN VISTA CTO4
POSTERIOR STRUT

MD MEDICAL DEVICE

SINGLE PATIENT - MULTIPLE USE

CAUTION








EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Pat. aspenmp.com/patents

Aspen Medical Products, LLC
6481 Oak Canyon, Irvine, CA 92618-5202 USA
aspenmp.com • +1-800-295-2776 • +1-949-681-0200

©2021 Aspen Medical Products, LLC. All Rights Reserved.
AC0130E revised 04/2021

ENGLISH
1. STRUT ASSEMBLY
NOTE: If MultiPost Collar has been applied, while maintaining proper spinal motion restriction protocols, remove the posterior panel and discard. Replace the posterior panel with the one provided and attach the posterior strut (see MultiPost instructions at www.aspenmp.com).
 A. Align three connection posts with three openings in MultiPost back panel and push into position.
 B. Push down causing strut support to lock into place. Ensure that all three connectors are engaged.
 C. Slide posterior strut into slot on lumbar belt posterior panel until it snaps into place.
2. CTO 4 POSTERIOR APPLICATION
NOTE: Maintain neutral alignment, and with proper support, log roll patient onto their side.
 A. Position posterior assembly and collar back panel on patient.
 B. To adjust the length push down on button. On backside, with other hand, pry back to fully release. Adjust strap up or down.
 C. If necessary, bend metal strap to account for kyphosis.
NOTE: Confirm adjustment buttons are secure following any adjustment.
 D. Log roll patient back onto their back.
 E. Tighten both collar straps as detailed in Vista MultiPost instructions.
 F. Snap lumbar straps into vest and tighten as needed to achieve a snug, comfortable fit.
NOTE: For Vista CTO 2 application, see Vista CTO 2 instructions www.aspenmp.com.

INTENDED USE
To provide motion restriction of the cervical spine.

INDICATIONS
Decompressive procedures, Fracture management, Postsurgical stabilization, Wiplash, Sprains/strains, Cervicogenic headache, Cervical disc syndrome, Musculoskeletal pain, Radiculopathy.

CONTRAINDICATIONS
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
Hand wash only at 30° C; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS
ABS, Nylon, HDPE, Polycarbonate, Delrin, Polyurethane, Cotton, Aluminium.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
Product may be safely disposed in accordance with local laws.

REPORTING
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
 (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
 (b) The product is altered;
 (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
 (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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(d) Le produit n'est pas utilisé conformément au mode d'emploi d'Aspen Medical Products, LLC.