
ASPEN HORIZON™ 639 LSO

- LUMBAR SACRAL ORTHOSIS
- LUMBAL-SAKRALE ORTHOSE
- ORTESIS LUMBOSACRA
- ORTHÈSE LOMBO-SACRÉE
- LUMBALE SACRALE ORTHÈSE
- BEĐERNÍ KRÍŽOVÁ ORTÉZA
- LUMBAL SAKRAL ORTÓZ
- ΟΡΘΩΣΗ ΟΡΘΩΣΗ ΙΕΡΟΥ ΟΣΤΟΥ
- LUMBAL KORSRYGGORTOSE
- ORTEZA ŁĘDŹIOWO-KRZYŻOWA
- LANNERANGAN SAKRALAINEN ORTOOSI
- LUMBAL SAKRAL ORTOS
- ОРТЕЗА ДЛЯ ПОЯСНИЧНО-КРЕСТЦОВОГО ОТДЕЛА ПОЗВОНОЧНИКА
- ÓRTESE SACRO LOMBAR
- 腰骶矫形器

1

2A

2B

2C

2D

2E

2F

3A

3B

4A

4B

4C

5A

5B

6

ENGLISH

1. SIZE ADJUSTMENT
To change the size of the Horizon, peel up the inside end of the side panel on each side. Adjust to proper fit and reattach.

Size	Circumference
1	24 – 34 in
2	32 – 40 in
3	38 – 46 in
4	44 – 50 in

2. PREPARATION
A. Place the pull tabs 2 inches (5 cm) from the plastic in the back of the brace. Pull the ends of the brace away from each other ensuring that it is fully extended.
B. Confirm that the lateral panels extend from posterior to anterior panels.
C. For smaller patients, later panels may need to be trimmed down to one panel. Once trimmed, tack exposed plastic end into the sleeve.
D. For shorter patients, the top extension panel may be taken off by removing the plastic anterior panel from the sleeve. Slide the aluminum bar to unlock. Once removed, reinsert the plastic anterior panel back into the sleeve.
E. For taller patients, the top extension panel, the aluminum bar may be bent to accommodate patient anatomy.
F. Align and attach left side panel to rigid anterior panel to accommodate patient anatomy.

3. APPLICATION
A. Apply the fully extended brace making sure the back panels is centered and low on the back. Apply the brace tightly around the abdomen, the right side overlapping the left.
B. When applied properly, both the upper and lower edge of the brace should be snug. If there is any gap, reapply the brace and angle the ends to insure a snug fit both the upper and lower edge.
NOTE: For certain patients, the ends must be steeply angled.

4. TIGHTENING
A. Slip your thumbs through the holes in the pull tabs. Pull away from your body.
B. The left pull tab controls lower compression and the right pull tab controls the upper compression. Adjust the compression as needed.
C. When tightened, secure the pull tabs to the hook-and-loop in front of the brace.

5. ADJUSTMENT
A. Adjust the upper and/or lower compression by loosening or tightening the pull tabs as needed.
B. When patient is seated, slightly loosening the upper portion may be beneficial.

6. MODIFICATION OPTIONS
This product may be substantially modified by an individual with expertise. For more information, please visit www.aspemp.com/education/product-modifications.

INTENDED USE
To provide motion restriction and stability of the lumbar-sacral spine.

INDICATIONS
Postoperative support, Deconditioned trunk musculature, Post laminectomy syndrome, Spinal Decompression procedures, Non-surgical Decompression procedures, Degenerative disc disease, Bulging or herniated disc, Non-specific low back pain, Chronic low back pain, Fracture management, Spondylolysis, Spondylolisthesis, Spinal stenosis, Ligament sprains/ Muscle strains, Spondylosis (Osteoarthritis/Degenerative Joint Disease), Facet syndrome, SI dysfunction, Sacroiliitis, Radiculopathy, Injections, Radiofrequency Ablation, Lumbar Multifidus Muscle Atrophy.

CONTRAINDICATIONS
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
Wash only at 30° C. Wash with mild detergent; Air dry. Do not use bleach or other chemicals.

MATERIALS
Nylon, Polyester, HDPE, PVC, Polyurethane, UHMWPE, Aluminum.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
This product may be safely disposed in accordance with local laws.

REPORTING
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts.

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
(a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
(b) The product is altered;
(c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
(d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

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