ASPEN SIERRA™ TLSO ADJUSTABLE

THORACIC LUMBAR SACRAL ORTHOSIS
LUMBAL-SAKRAL-THORAKALE ORTHESE

ÓRTEZIS TORÁCICA LUMBOSACRA

ORTOSI TORACO-LOMBO-SACRALE

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LUMBALE SACRALE THORACALE ORTHESE

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









ÓRTEZA TORÁCICA SACRO LOMBAR

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ENGLISH

1. SIZE ADJUSTMENT
A. To adjust the size of the Aspen Sierra TLSO, peel up the inside end of the side panel on each side. Adjust to achieve a proper fit and reattach.
B. Size ranges are only guidelines. Adjust as needed to achieve a snug fit without too much overlap.

Size	Circumference
1	28 – 37 in 71 – 94 cm
2	36 – 43 in 91 – 109 cm
3	42 – 48 in 107 – 122 cm
4	47 – 54 in 119 – 137 cm

2. PREPARATION
Place the pull tabs 2 inches (5 cm) from the plastic in the back of the Aspen Sierra TLSO brace. Pull the ends of the brace away from each other ensuring that it is fully extended.

3. LATERAL PANEL SIZING
Add or remove panels for appropriate sizing.

4. PLACEMENT & ADJUSTMENT
A. Center the back panel as low as possible. Slide abdominal panel and sternal pad sideways as needed to ensure brace and panel are at patient's midline.
B. To determine sternal pad height, pull button back and slide aluminum strut so that the button locks into the appropriate hole.

5. TIGHTENING
Slip your thumbs through the holes in the pull tabs. Pull away from your body. The pull tab controls lower compression and the right pull tab controls the upper compression. Adjust the compression as needed. When tightened, secure the pull tabs to the loop in front of the brace.

6. SHOULDER STRAPS
A. Attach control straps to back panel. Straps can be positioned to fit under the arms or over the shoulders.
B. Control straps and sternal pad should fit securely to patient's anatomy. Pull straps to tighten as needed.
C. Torso control strap setting.
D. Over the shoulder setting.

7. MODIFICATION OPTIONS
This product may be substantially modified by an individual with expertise. For more information, please visit www.aspenmp.com/education/product-modifications.

INTENDED USE
To provide motion restriction and stability of the Thoracic-lumbar-sacral spine.

INDICATIONS
Postoperative support, Degenerative disc disease, Bulging or herniated disc, Fracture management, Spondylolysis, Spondylolisthesis, Spinal stenosis, Kyphosis, Spondylolysis (Osteoarthritis/Degenerative Joint Disease).

CONTRAINDICATIONS
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
Hand wash only at 30° C; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS
Nylon, Polyester, HDPE, PVC, Polyurethane, UHMWPE, Aluminum.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
Product may be safely disposed in accordance with local laws.

REPORTING
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OR DEFECT IN THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALSO OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

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Dans les cas suivants, Aspen Medical Products, LLC n'aura aucune obligation au titre de cette garantie limitée :
(a) Le produit n'a pas été acheté auprès d'Aspen Medical Products, LLC ou via ses canaux de distribution autorisés ;
(c) Toute pièce n'est fournie par Aspen Medical Products, LLC à l'exception de celle dans le produit ; ou
(d) Le produit n'est pas utilisé conformément au mode d'emploi d'Aspen Medical Products, LLC.