








**OSTEOARTHRITIS KNEE BRACE**  
**OSTEOARTHRITIS-KNIEORTHESE**  
**ÓRTEZIS DE RODILLA PARA OSTEOARTRITIS**  
**GINOCCHIERA PER L'OSTEOARTRITE**  
**ORTHÈSE POUR ARTHROSE DU GENOU**  
**KNIEBRACE TEGEN ARTROSE**  
**KOLENNÍ ORTÉZA PRO OSTEOARTRÓZU**  
**KNÆBØJLE TIL SLIDGIGT**  
**ΝΑΡΘΗΚΑΣ ΓΟΝΑΤΟΥ ΓΙΑ ΟΣΤΕΟΑΡΘΡΙΤΙΔΑ**  
**KNESTØTTE FOR OSTEOARTRITT**  
**ORTEZA KOLANOWA DLA PACJENTÓW Z CHOROBY ZWYRONNIENIOWĄ STAWÓW**  
**POLVEN NIVELRIKKO-ORTOOSI**  
**KNÄSTÖD FÖR ARTROS**  
**Коленный ортез при остеоартрозе**  
**JOELHOIRA PARA OSTEOARTRITE**  
**骨关节炎膝关节支架**

Aspen Medical Products, LLC  
 6481 Oak Canyon, Irvine, CA 92618-5202 USA  
 aspenmp.com • +1-800-295-2776 • +1-949-681-0200  
 ©2021 Aspen Medical Products, LLC. All Rights Reserved.  
 LB0051B revised 12/2021

**ENGLISH**  
**1. FITTING THE BRACE**  
 A. Before fitting the brace, ensure all brace fasteners are open and set the horizontal position of the hinge to 0°.  
 B. While sitting, bend the knee joint slightly (approximately 30°).  
 C. Position the brace on the side of the leg so that the center of the hinge lines up with the center of the patella.  
 D. Close fasteners starting with the upper calf tab #1. Adjust and tighten the strap so that it is comfortable for the patient and prevents the brace from migrating. Follow the same process with lower thigh tab #2, then lower calf tab #3 and then upper thigh tab #4. Pull the calf strap #5 through the buckle and secure tightly. Repeat the process for the thigh strap #6.  
 E. Use the adjustment tool to set the appropriate degree of unloading. Clockwise decreases unloading and counterclockwise increases unloading. While the patient is sitting, turn the adjustment screw 1/8 to a 1/4 turn at a time and evaluate patient before adjusting further.

**2. DOFFING THE BRACE**  
 To doff the brace, loosen and remove the thigh #1 and calf strap #4. Remove the tabs and doff the brace.  
**3. ADJUSTING THE RANGE OF MOTION**  
 If required, flexion / extension limits can be set. Remove the condyle pad by sliding the pad out of the wrap material. Remove the screw for the limit stop. Remove the limit stop from the hinge and insert the required angular degree. Secure the inserted limit stop with the screw and reinsert the pad into the condyle pocket.  
**Limits available are:**  
 -Extension: 0°, 5°, 15°, 20° and 25°  
 -Flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°

**INTENDED USE**  
 To provide medial or lateral forces on the knee joint.  
**INDICATIONS**  
 Osteoarthritis, medial or lateral compartment of the knee.  
**CONTRAINDICATIONS**  
 Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.  
**PRECAUTIONS**  
 READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

**CARE INSTRUCTIONS**  
 Hand wash only at 30° C; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.  
**MATERIALS**  
 Spandex fiber, Aluminum, Stainless steel, Nylon, Carbon steel, Polyurethane, Polyester, EVA, Polypropylene.  
**NOT MADE WITH NATURAL RUBBER LATEX.**  
**DISPOSAL**  
 Product may be safely disposed in accordance with local laws.

**REPORTING**  
 Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).  
**LIMITED WARRANTY**  
 Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).  
 Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:  
 (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution.  
 (b) The product is altered;  
 (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or  
 (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.  
 THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.  
 A listing of authorized distributors in your area is available upon request.

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(c) Cualquier repuesto no suministrado por Aspen Medical Products, LLC fue insertado en el producto.
(d) El producto no se utiliza de acuerdo con las Instrucciones de Uso de Aspen Medical Products, LLC.