

Aspen Vista[®] MultiPost Therapy Collar

CODE L0180 APPROVED

The Aspen Vista MultiPost Therapy Collar is designed to address the symptoms associated with cervical conditions such as Forward Head Carriage ("Text Neck"). Unnatural head positions can shift the cervical spine out of alignment, leading to radiculopathy, neuropathy and cervicogenic pain. If left untreated, the unbalanced loading of the cervical spine can lead to permanent deterioration of discs, vertebrae, ligaments and muscles. This device actively influences posture and a more natural alignment to aid in the treatment of these disorders.

Active Influence from the Comfort of Home



Aspen Active Influence Braces are Designed to Promote a More Natural Alignment




Aspen's Active Influence line of braces are specifically designed to help patients suffering from chronic postural conditions. These innovative braces provide postural support with the goal of reducing symptoms¹ caused by an unbalanced body position. Aspen Active Influence braces also work to encourage symmetrical spinal loading, prompting a more natural center of gravity and a restoration of normal proprioception.

1. Fabio Zaina, Martina Poggio, Sabrina Donzelli, and Stefano Negrini. "Can Bracing Help Adults with Chronic Back Pain and Scoliosis? Short-Term Results from a Pilot Study." *Prosthetics and Orthotics International*. 2018.

The Aspen Advantage Symbol
An Innovative Feature Available only from Aspen

Treating Forward Head Carriage 

The collar applies gentle pressure to the chin in the front and to the mid cervical region in the back to counteract the unnatural forces caused by Forward Head Carriage. Worn in the privacy of one's home, this unique device is designed with the goal to help relieve cervicogenic headaches, radiculopathy and reduce pain.



One Size Adjustable 

Patented sizing system offers six height adjustment settings to easily accommodate a wide range of anatomies, allowing for easy adjustments on the patient.



One-Time Adjustment 

Once the Locking Strap on the left side is adjusted to the correct position, the patient can easily achieve the proper tightness and positioning when the collar is reapplied at home.

Integrated Air Pump 

Easily accessible, the air pump and pressure release valve allow patients to quickly achieve the level of compression that works best for them.




NEWMATIX™
COMPONENT
THERAPY
PACK



Active Influence Compression 

This device provides therapeutic support to actively influence a more natural curve in the neck. Activating the NEWMatix™ System bladder can aid in healing and help reduce symptoms by taking stress off the vertebra, discs and nerves.



Targeted Therapy 

The Aspen Therapy Pack delivers hot/cold therapy to reduce inflammation and spasm. Moving its location inside the back panel focuses the targeted compression where desired.

Aspen Vista® MultiPost Therapy Collar - Code L0180 Approved

984250

Aspen Vista® MultiPost Therapy Collar Replacement Pad Set

984252

Aspen Vista® MultiPost Therapy Collar Replacement Chin Pad

984255

Aspen Vista® MultiPost Therapy Collar Replacement Therapy Pack

984257

