

Aspen Vista[®] Cervical Collar

CODE L0172 APPROVED

The fully adjustable Aspen Vista Cervical Collar helps improve patient care while saving time, money and storage space. Built on the design of the trusted Aspen Cervical Collar, the Aspen Vista Collar provides motion restriction¹ while helping to protect against skin breakdown. With its patented height adjustment technology, the right size is always on hand.

Award-Winning Cervical Care

The unique and innovative design of the Aspen Vista Collar earned both the prestigious *SPINE TECHNOLOGY AWARD* for Best New Technology in Cervical Care and the TechAmerica *TERMAN AWARD* for the best Medical Device. Aspen is proud to have won these awards and to be the cervical orthosis of choice in the majority of Trauma Centers across North America, working together to provide superior patient care and improved outcomes.



TechAmerica




Achieve a Controlled Environment for Healing



The Aspen line of Acute Restriction collars is designed to achieve a controlled environment for patients healing from trauma or surgery. These collars provide optimum cervical motion restriction², minimize incidences of skin breakdown³ and increase patient compliance.

- 1, 2. Ivancic, Paul C. "Do Cervical Collars and Cervicothoracic Orthoses Effectively Stabilize the Injured Cervical Spine? A Biomechanical Investigation." *Spine*, vol. 38, no. 13, 2013, pp. 767-74.
3. Powers, Jan, et al. "The Incidence of Skin Breakdown Associated With Use of Cervical Collars." *Journal of Trauma Nursing*, vol. 13, no. 4, 2006, pp. 198-200.

The Aspen Advantage Symbol
An Innovative Feature Available only from Aspen



One Size Adjustable
Patented sizing system offers six height adjustment settings to easily accommodate a wide range of anatomies, allowing for easy adjustments on the patient.



Dial Height Adjustment [Ⓐ]
Simply pull and turn the Height Adjustment Dial to select one of six height settings.



Aspen Vista® TX Collar
CODE L0174 APPROVED

Encompasses all of the features and benefits of the popular Aspen Vista Cervical Collar with an easily detachable thoracic extension, providing added coverage in the thoracic region.



Extra Large Access Window [Ⓐ]
The extra large tracheal aperture provides functional access for airway management and other essential procedures.

Aspen Vista® Collar	984000
Aspen Vista® Collar Set with Replacement Pads	984002
Aspen Vista® Collar Replacement Pads	984020
Aspen Vista® TX Collar	984100
Aspen Vista® TX Collar Set with Replacement Pads	984102
Aspen Vista® TX Collar Replacement Pads	984120
Standard Back Panel	984006
Large Back Panel	984008
Large Back Panel Replacement Pad	984028