



# WALK A MILE IN THEIR SHOES. (OR 2000)

Having undergone 2 million cycles of testing, the M2 AFO™ has proven to be durable for extended use to facilitate a normal gait. The lightweight and comfortable design offers ideal support with optimal flex for toeoff.





# M2 AFO™



## Carbon Fiber

Creates dynamic forces for comfortable walking.




## Diamond Stitched Fabric

Provides comfort for all-day wear.



## Trimmable Footplate

Can be customized to fit most patient anatomies and shoes.



## Medium Flex

Fits most patient needs for enhanced performance.



	POSTERIOR M2 AFO - Code L1951		ANTERIOR M2 AFO - Code L1932	
SIZE	LENGTH	HEIGHT	LENGTH	HEIGHT
S	10.04 in (255 mm)	13.78 in (350 mm)	10.04 in (255 mm)	15.35 in (390 mm)
M	10.63 in (270 mm)	14.17 in (360 mm)	10.63 in (270 mm)	15.75 in (400 mm)
L	11.42 in (290 mm)	14.57 in (370 mm)	11.42 in (290 mm)	16.14 in (410 mm)
XL	11.61 in (295 mm)	14.96 in (380 mm)		

	RIGHT		LEFT	
POSTERIOR	Small - Code L1951	640112	Small - Code L1951	640212
	Medium - Code L1951	640113	Medium - Code L1951	640213
	Large - Code L1951	640114	Large - Code L1951	640214
	XLarge - Code L1951	640115	XLarge - Code L1951	640215
ANTERIOR	Small - Code L1932	640102	Small - Code L1932	640202
	Medium - Code L1932	640103	Medium - Code L1932	640203
	Large - Code L1932	640104	Large - Code L1932	640204
ACCESSORIES	Anterior Replacement Pads, Right		640101	
	Anterior Replacement Pads, Left		640201	
	Posterior Replacement Pads, Universal		640111	