

Aspen Hinged Knee

Maximum Knee Support In A Low-Profile Design

XS	648001	L	648004
S	648002	XL	648005
M	648003	XXL	648006

The Aspen Hinged Knee sized brace features a premium AirX material that provides increased compression and support with a high degree of breathability and moisture wicking. The comfort provided by the AirX material also allows for extended wear and improves patient compliance.

Aspen Hinged Knee Universal




Provides adjustability to fit most patients. Made from a lightweight neoprene material, this option allows for optimal compression and comfort.

Aspen Hinged Knee Universal
Code L1820 Approved

648011



The perfect blend of comfort and control.



Comfort

Advanced AirX material provides ideal compression with a high degree of breathability and moisture wicking.



Hinge

Polycentric hinge for peak performance.



Pull Tabs

Molded pull tabs for easy adjustment.



Condyle Pads

Integrated condyle pads for added comfort.



Aspen Knee Suspension Wrap

The Aspen Knee Suspension Wrap can be used in conjunction with the Aspen Hinged Knee for patients needing additional support.

Accessory

Aspen Knee Suspension Wrap
Code L2397 Approved

651111

Aspen Hinged Knee XS - Code L1820 Approved - 11-15 in (27.94-38.1 cm)	648001
Aspen Hinged Knee S - Code L1820 Approved - 15-18 in (38.1-45.72 cm)	648002
Aspen Hinged Knee M - Code L1820 Approved - 18-21 in (45.72-53.34 cm)	648003
Aspen Hinged Knee L - Code L1820 Approved - 21-24 in (53.34-60.96 cm)	648004
Aspen Hinged Knee XL - Code L1820 Approved - 24-27 in (60.96-68.58 cm)	648005
Aspen Hinged Knee XXL - Code L1820 Approved - 27-35 in (68.58-88.9 cm)	648006
Aspen Hinged Knee Universal - Code L1820 Approved	648011