

# WE'VE CREATED A BRACE THAT'S PERSONALIZED TO YOU.

A streamlined design, customizability, and comfort create an easy, all-day wearing experience – because you deserve an individualized approach to healing.



Aspen  
MEDICAL PRODUCTS



## SIERRA™ HYPEREXTENSION TLSO

The All Day Solution to Hyperextension.

Aspen  
MEDICAL PRODUCTS

© 2021 Aspen Medical Products, LLC  
800-295-2776 • 949-681-0200 • aspenmp.com  
All rights reserved. Pat. aspenmp.com/patents.  
LIT5170 12/21



# DESIGNED WITH YOUR COMFORT IN MIND.




## TIGHTENING CONTROL

Provides the ideal amount of tightening after putting on the brace.



## PIVOTING PELVIC BAR

Prevents the brace from riding upward when you move from standing to sitting.



## BREATHABLE MATERIAL

For a comfortable healing experience.



SCAN HERE WITH YOUR  
PHONE CAMERA TO WATCH  
A VIDEO ABOUT HOW TO  
APPLY, REMOVE, AND CARE  
FOR YOUR BRACE.

