

# Traverse Air Walker - Short

The Aspen Traverse Air Walker is a circumferential walking boot that utilizes dual air pumps for targeted, comfortable compression. Walker boots are designed to help prevent further damage and promote quicker healing of the injured site. The walker boot helps the user to become mobile again by protecting the site and slowly allowing weight to be put onto the injured foot/ankle.

- Comfortable and lightweight
- Easy to apply and remove
- Air pump pressure can be controlled from either side of the ankle
- Rocker base promotes a natural and stable gait
- Shock absorbing heel design



**CODE L4360/L4361 APPROVED**

Size	Men's Shoe Size	Women's Shoe Size
XS	Up to 3.5	Up to 4.5
S	3.5 - 6.5	4.5 - 7.5
M	6.5 - 9	7.5 - 10
L	9 - 11	10 - 12
XL	11 - 14.5	12 - 16

Size					Accessories				
XS	S	M	L	XL	Achilles Wedge Kit	RIGHT	JWW/RR	LEFT	JWW/RL
SKU									
653001	653002	653003	653004	653005	Extended Strap Kit				ASW/EXS

Traverse Replacement Liner — Short					Traverse Hygiene Cover				
XS	S	M	L	XL	XS	S	M	L	XL
653301	653302	653303	653304	653305	653401	653402	653403	653404	653405

# Traverse Air Walker - Tall

The Aspen Traverse Air Walker is a circumferential walking boot that utilizes dual air pumps for targeted, comfortable compression. Walker boots are designed to help prevent further damage and promote quicker healing of the injured site. The walker boot helps the user to become mobile again by protecting the site and slowly allowing weight to be put onto the injured foot/ankle.

- Comfortable and lightweight
- Easy to apply and remove
- Air pump pressure can be controlled from either side of the ankle
- Rocker base promotes a natural and stable gait
- Shock absorbing heel design



## CODE L4360/L4361 APPROVED

Size	Men's Shoe Size	Women's Shoe Size
XS	Up to 3.5	Up to 4.5
S	3.5 - 6.5	4.5 - 7.5
M	6.5 - 9	7.5 - 10
L	9 - 11	10 - 12
XL	11 - 14.5	12 - 16

Size					Accessories				
XS	S	M	L	XL	Achilles Wedge Kit	RIGHT	JWW/RR	LEFT	JWW/RL
SKU									
653011	653012	653013	653014	653015	Extended Strap Kit				ASW/EXS

Traverse Replacement Liner — Tall					Traverse Hygiene Cover				
XS	S	M	L	XL	XS	S	M	L	XL
653306	653307	653308	653309	653310	653401	653402	653403	653404	653405