

Elite™ TLSO

CODE L0456 / L0457 APPROVED

The Elite TLSO provides clinicians with a light weight, low profile option to treat patients with symptoms related to osteoporosis. The moldable, adjustable, rigid posterior strut combined with an inelastic knitted belt helps to address kyphosis which could lead to better balance and more daily activities. The patented independent compression mechanism allows patients to target and adjust the compression for optimum comfort and pain relief.





The Aspen Advantage Symbol
An Innovative Feature Available Only From Aspen

Height Adjusts for Every Body
Telescoping shoulder straps provide 6 inches of height adjustment to fit varying patient anatomies.

Adjustable Strut Adjustable and customizable strut height to fit all patient sizes.

UltraKnit™ Technology

Knitted technology integrates both elastic and inelastic components.



Independent Tightening Mechanism

Independent compression for ideal pain relief.

Customizable straps

Customizable strap length for optimum support and comfortable wear.



Removable strut

Removable strut allows for use of lumbar belt to be used independently.

← FITS WAIST CIRCUMFERENCE RANGING FROM 25.5-59.05 IN (65-150 CM) AND UP TO 68.5 IN (175CM) WITH TLSO EXTENSION PANEL →

Elite™ TLSO – XS, Code L0456 / LO457 Approved - 25.5-29.5 in (65-75 cm)	591601
Elite™ TLSO – S, Code L0456 / LO457 Approved - 29.5-33.5 in (75-85 cm)	591602
Elite™ TLSO – M, Code L0456 / LO457 Approved - 33.5-37.25 in (85-95 cm)	591603
Elite™ TLSO – L, Code L0456 / LO457 Approved - 37.25-41.25 in (95-105 cm)	591604
Elite™ TLSO – XL, Code L0456 / LO457 Approved - 41.25-45.25 in (105-115 cm)	591605
Elite™ TLSO – XXL, Code L0456 / LO457 Approved - 45.25-49.25 in (115-125 cm)	591606
Elite™ TLSO – XXXL, Code L0456 / LO457 Approved - 49.25-59.05 in (125-150 cm)	591607
Elite™ TLSO – Extension Panel - Added 9.5 in (25 cm)	590555