

# VRTX TLSO & LSO

The VRTX TLSO has integrated years innovations into an intuitive and specialized design for efficient support within the thoracolumbar spine, providing patients with an innovative brace that delivers appropriate motion restriction throughout the healing process. The VRTX TLSO's simplified two piece bivalved design allows for clinicians to apply the brace in a supine position reducing the amount of patient log rolling. Improved under arms straps provide optimal stabilization and comfort. This one size fits most brace utilizes built in side panels that extend to fit a wide variety of patient anatomies. This key feature assists in reducing part inventory for clinics. The VRTX TLSO can be easily stepped down to an LSO or stepped up to a CTLSO as part of the interactivity found in the VRTX family of products.




## VRTX TLSO

Ideal motion restriction for acute thoracolumbar spine trauma and post-operative patients.  
Code L0464 approved



## VRTX LSO

As the patient improves and treatment plan progresses, the VRTX TLSO has the ability to step down to an LSO to provided the needed level of support that fits into your practices' plan of care.  
Code L0637 approved




# Direct customizable compression and increased trunk stability.



**Optimal Motion Restriction**  
Underarm straps for optimal patient motion stabilization.



**Two Piece Design**  
Simplified two-piece bivalved design for easier patient re-application.



**Modularity In Treatment**  
The VRTX TLSO can be stepped down to an LSO or CTLSO based on patient's need.



**Versatility In Sizing**  
Fit a wide variety of patients with the VRTX TLSO specially designed anterior panel.



**Limit Logrolling**  
Specialized anterior tightening and adjustments allows for more efficient clinician application.

## TLSO

VRTX TLSO - Standard	995610	VRTX TLSO - XL	985630
VRTX TLSO - XS	995620		

## LSO

VRTX LSO - Standard	995410	VRTX LSO - XS	995420
VRTX LSO - XL	995430		

## ACCESSORIES

VRTX Lower Spine Replacement Pads - Standard	995613	VRTX Lower Spine Replacement Pads - XL	995633
VRTX Lower Spine Replacement Pads - XS	995623	VRTX Shoulder Strap Kit	995002