

Symmetrical Back Panel

The Vista II collar features an innovative symmetrical back panel that integrates with the VRTX CTLSO for a cervical to sacrum solution.




Vista II Cervical Collar

Incorporates the trusted design of the original Vista Cervical Collar, including a one size adjustable fit and dial height adjustment. Code L0172 approved

Vista II Cervical Collar & Vista II CTLSO

The Vista II Cervical Collar integrates the proven technology of the original Vista Cervical Collar with the reversible back panel of the VRTX Collar, ensuring precise and consistent application. This not only optimizes application times but also proves invaluable in critical situations, promoting efficiency and reliability in clinical settings.

This comprehensive off the-shelf cervical to sacrum solution enhances the overall efficiency of hospital operations.



Symmetrical Back Panel

The Vista II collar features an innovative symmetrical back panel that streamlines fitting and reduces the risk of incorrect patient application. The symmetrical posterior panel serves as a visual cue, ensuring that the collar is applied accurately and securely.




Connecting Straps

Adjustable straps easily connects the Vista II collar with the Vista II CTLSO upgrade kit for multilevel spinal support.



Cervical and Sternal Adjustability

Easy multi step cervical and sternal adjustability allows clinician to provide the care and comfort needed for the patient.



Vista II CTLSO Upgrade Kit

For multilevel spinal support, step up the Vista II easily with the Vista II CTLSO upgrade kit.



VRTX LSO/TLSO

Mate the Vista II Collar with the VRTX LSO/TLSO and the Vista II CTLSO upgrade kit for a comprehensive cervical to sacrum solution.

VISTA II COLLAR

Vista II Collar	984300
Vista II Collar Set with Replacement Pads	984302
Vista II Collar Replacement Pads	984320

VISTA II CTLSO

Vista II CTLSO Upgrade Kit - Standard/XL	984650
Vista II CTLSO Upgrade Kit - XS	984660

VISTA II CTLSO STANDALONE

Vista II CTLSO - Standard	984680
Vista II CTLSO - XS	984690
Vista II CTLSO - XL	984700