

## Aspen Vista<sup>®</sup> CTO and Aspen Vista<sup>®</sup> CTO4

Designed to provide motion restriction throughout the cervico-thoracic spine, the Aspen Vista CTO and Aspen Vista CTO4 provide significant motion control in all three planes of motion—flexion/extension, lateral bending and axial rotation.<sup>1</sup> These one size adjustable braces feature multiple adjustment points to help ensure a comfortable fit and the ability to step down to support patient recovery.

### Award-Winning Cervical Care

The unique and innovative design of the Aspen Vista CTO earned the prestigious *SPINE TECHNOLOGY AWARD* for Best New Technology in Cervical Care. Aspen is proud to have won this award and to be the cervical orthosis of choice in the majority of Trauma Centers across North America, working together to provide superior patient care and improved outcomes.



## Achieve a Controlled Environment for Healing



The Aspen line of Acute Restriction collars is designed to achieve a controlled environment for patients healing from trauma or surgery. These collars provide optimum cervical motion restriction<sup>2</sup>, minimize incidences of skin breakdown<sup>3</sup> and increase patient compliance.

- 1, 2. Ivancic, Paul C. "Do Cervical Collars and Cervicothoracic Orthoses Effectively Stabilize the Injured Cervical Spine? A Biomechanical Investigation." *Spine*, vol. 38, no. 13, 2013, pp. 767–74.
3. Powers, Jan, et al. "The Incidence of Skin Breakdown Associated With Use of Cervical Collars." *Journal of Trauma Nursing*, vol. 13, no. 4, 2006, pp. 198–200.


**The Aspen Advantage Symbol**  
An Innovative Feature Available only from Aspen



**Padded for Comfort**  
Fully padded at all patient contact points to help maximize comfort, encouraging improved patient compliance.



**Locking Sternal Adjustment** <sup>Ⓐ</sup>  
Angle of support strut is easily adjusted to fit a patient's anatomy then locks into place.



**Flexible Thoracic Adjustment** <sup>Ⓐ</sup>  
One size fits essentially all patient anatomies and unlocks to adjust during original fit and repositioning.



ASPEN VISTA CTO4



ASPEN VISTA CTO



ASPEN VISTA MULTIPOST COLLAR



**The Aspen Vista CTO4** <sup>Ⓐ</sup>  
For enhanced motion restriction, the Aspen Vista CTO4 provides increased motion control in all three planes of motion. The adjustable, low-profile back strut is designed for maximum patient comfort, even in the supine position.

**Easily Steps Down to Support Progression of Care** <sup>Ⓐ</sup>

As a patient's condition improves, the Aspen Vista CTO4 can be stepped down to the Aspen Vista CTO then down again to the Aspen Vista MultiPost Collar, all without the use of tools.

Aspen Vista® CTO4 - 4 Post CTO	984550
Aspen Vista® CTO - 2 Post CTO	984510
Aspen Vista® CTO4 Upgrade Kit - Posterior	984554
Aspen Vista® CTO Upgrade Kit - Anterior	984504
Aspen Vista® CTO4 / Vista CTO Replacement Pads	984520