

Peak Scoliosis Bracing System

The Peak Scoliosis Bracing System helps to provide relief of the symptoms associated with Adult Scoliosis. This Active Influence system uses a combination of compression and support to promote better alignment, improve posture and restore a more natural center of gravity.

Designed for Unique Patient Needs

Multiple configuration options give providers the ability to customize a patient's treatment. With this modular design, providers are able to assemble the right brace for each patient.




MULTIPLE OPTIONS FROM ONE BRACE

DUAL STRUT
CONFIGURATION KIT



Address pain due to adult scoliosis for improved activities of daily living.



SlickTrack™ Tightening System
Features an independent upper and lower tightening mechanism, providing localized compression where needed.

Convenient Closure Handles
Enhances ease of use when applying or removing the brace.



Enhanced Pulley System
Reduces friction and makes it easy to control the level of compression regardless of patient's strength.

Adjustable Tension Straps

To optimize the angle from which tension is applied, the posterior tension strap can be placed anywhere along the attachment strips based on patient needs.



Foam for Comfort

Reliable Structural Integrity
Integrates vertical stays that provide stability and can be removed on smaller sizes for a low profile, comfortable fit.

Universal Strut Assembly

Can be adjusted or modified to accommodate patient's unique needs, then placed anywhere on the belt to optimize support.

← FITS WAIST CIRCUMFERENCE RANGING FROM 24-45 IN (61-114 CM) AND UP TO 55 IN (139 CM) WITH ASPEN EXTENSION PANEL →



- INCLUDES:**
- Peak Scoliosis Belt
 - Universal Strut Assembly
 - Lateral Panel with Iliac Pad
 - Chest Strut & Chest Pad
 - Carry Bag
 - Instructions for Use
 - Patient Handbook

One Size Adjustable - 24-45 in (61-114 cm)		993911
Small - 24-30 in (61-76 cm)	993901	Dual Strut Configuration Kit 993905
Medium - 29-35 in (74-89 cm)	993902	Peak Short Trochanter Strut 993916
Large - 34-40 in (86-102 cm)	993903	Replacement Sleeve Set 993910
X-Large - 39-45 in (99-114 cm)	993904	Aspen Extension Panel - 10 in (25 cm) 992505