

Aspen Contour™ Bracing System

The Aspen Contour Bracing System provides multiple orthotic configurations in one system, offering an effective treatment solution to support various clinical indications. This product line delivers a comfortable fit, effective compression and the rigidity required to create a controlled healing environment for pre-, post- and non-operative patients.

Custom Build the Right Brace

Practitioners can custom build the right brace to meet each patient's clinical needs and unique anatomy through various configuration options. To customize the brace further, various components can be removed, added or significantly modified.



Aspen Acute Restriction Braces Provide a Controlled Environment for Healing



Aspen Acute Restriction braces provide significant motion restriction by acting as a kinematic restrictor,¹ for patients needing support for healing, and can be safely and quickly applied. Gross motion restriction is often desired in treating injured patients recovering from trauma or surgery.² The goal is to support healing,³ limit recovery time and save money by expediting patient discharge times, accelerating turnaround times and minimizing health risks.

1. Cholewicki J, Lee AS, Reeves NP, Morrissette DC. "Comparison of trunk stiffness provided by different design characteristics of lumbosacral orthoses." *Clin Biomech (Bristol, Avon)*. 2010 Feb; 25(2): 110–114., doi: 10.1016/j.clinbiomech.2009.10.010.

2. Cholewicki, Jacek, et al. "Comparison of Motion Restriction and Trunk Stiffness Provided by Three Thoracolumbosacral Orthoses (TLSOs)," *Journal of Spinal Disorders & Techniques*, vol. 16, no. 5, Oct. 2003, pp. 461–468.

3. Cholewicki, Jacek. "The Effects of Lumbosacral Orthoses on Spine Stability: What Changes in EMG Can Be Expected?" *Journal of Orthopaedic Research*, vol. 22, no. 5, 2004, pp. 1150–1155., doi:10.1016/j.jorthres.2004.01.009.

The Aspen Advantage Symbol

An Innovative Feature Available only from Aspen

Custom Build the Right Brace

The Contour Bracing System can be easily stepped up or down to accommodate varying patient conditions and to further support the continuum of care.

Easy Access for Modification

All structural components are easily accessible and can be heat molded, bent or trimmed to accommodate individual patient conditions.



CONTOUR TLSO CHEST PANEL




RIGID ANTERIOR PANEL (RAP)



ASPEN QUIKDRAW™ PRO

CONTOUR LSO BACK PANEL



CONTOUR TLSO SHOULDER STRAPS




Custom Support

When tightened, the firm but flexible QuikDraw conforms to the patient's lordosis to provide true structural support.

Effective Independent Compression

Aspen's patented 4:1 mechanically advantaged tightening system independently tightens top and bottom; allowing patients to localize compression where they need it.



Aspen Contour™ TLSO SPK Complete - Code L0464 Approved							
Medium - 31-37 in (79-94 cm)	994351	Large - 36-42 in (91-107 cm)	994352	X-Large - 41-47 in (104-119 cm)	994353	XX-Large - 46-57 in (117-145 cm)	994354
TLSO SPK Upgrades - Small	994455	TLSO SPK Upgrades - Medium	994460	TLSO SPK Upgrades - Large	994465	TLSO SPK Replacement Pads	994470
Contour TLSO Upgrades - Small Chest Panel			994255	Contour TLSO Upgrades - Medium Chest Panel			994260
Contour TLSO Upgrades - Large Chest Panel			994265	Contour TLSO Upgrades Short Shoulder Strap	994151	Contour TLSO Upgrades Tall Shoulder Strap	994152
Aspen Contour™ LSO Complete - Code L0637 Approved							
X-Small - 21-27 in (53-69 cm)	990050	Medium - 31-37 in (79-94 cm)	990052	X-Large - 41-47 in (104-119 cm)	990054	12" Contour Back Panel	992212
Small - 26-32 in (66-81 cm)	990051	Large - 36-42 in (91-107 cm)	990053	XX-Large - 46-57 in (117-145 cm)	990055	15" Contour Back Panel	992215
12" Back Panel Replacement Pad			992122	15" Back Panel Replacement Pad			992125
Aspen QuikDraw™ RAP - Code L0627 Approved							
Black X-S - 21-27 in (53-69 cm)	995010	Black Large - 36-42 in (91-107 cm)	995040	White X-S - 21-27 in (53-69 cm)	995011	White Large - 36-42 in (91-107 cm)	995041
Black Small - 26-32 in (66-81 cm)	995020	Black X-L - 41-47 in (104-119 cm)	995050	White Small - 26-32 in (66-81 cm)	995021	White X-L - 41-47 in (104-119 cm)	995051
Black Medium - 31-37 in (79-94 cm)	995030	Black XX-L - 46-57 in (117-145 cm)	995060	White Medium - 31-37 in (79-94 cm)	995031	White XX-L - 46-57 in (117-145 cm)	995061
Aspen QuikDraw™							
Black X-S - 21-27 in (53-69 cm)	993081	Black Large - 36-42 in (91-107 cm)	993031	White X-S - 21-27 in (53-69 cm)	993080	White Large - 36-42 in (91-107 cm)	993030
Black Small - 26-32 in (66-81 cm)	993011	Black X-L - 41-47 in (104-119 cm)	993041	White Small - 26-32 in (66-81 cm)	993010	White X-L - 41-47 in (104-119 cm)	993040
Black Medium - 31-37 in (79-94 cm)	993021	Black XX-L - 46-57 in (117-145 cm)	993091	White Medium - 31-37 in (79-94 cm)	993020	White XX-L - 46-57 in (117-145 cm)	993090
Black Regular RAP	993200	Black Large RAP	993201	White Regular RAP	993100	White Large RAP	993101
Regular RAP Replacement pad	993110	Large RAP Replacement Pad	993111	Black Extension Panel	993061	White Extension Panel	993060
Cold Therapy Pack							993005