



Aspen[®]
LSO

High Performance - Low Profile

The Aspen[®] LSO Low Profile provides immediate pain relief and comfortable support. When tightened using the LeverLatch[™] mechanism, the brace firmly and evenly compresses the torso. This enables patients to regain a significant degree of functional mobility and, in many cases, return to routine activities of daily living.

The unique, three-piece design, featuring symmetrical side panels, allows the practitioner to configure the orthosis so it can be tightened using either hand.

When coupled with cotton-lined, breathable replacement pads and CushionFlex Tab[™] edges, the Aspen[®] LSO comfortably "hugs" the torso and can be easily worn under a jacket or sweater.

Pivoting LeverLatch[™] Tightener

Allows the brace to self-adjust to various physiques for a better fit. Offers a powerful mechanical advantage to apply even compression for immediate pain relief.

Cylindrical Support Structure

Encircles the torso with a structurally rigid polyethylene cylinder for strong support and better immobilization, yet is easy to trim or shape.

Comfortable CushionFlex Tabs[™]

Softens the edges of the orthosis and conform to the curves of the patient for a comfortable fit and better patient compliance.



Innovative Patient Care

Effective Versatile Bracing

By choosing to use an Aspen[®] LSO early in the care of patients with lumbosacral problems, practitioners can make an immediate and positive impact.

Simply form the back panel, adjust the sides, tighten for support, and trim for comfort.

Whether used as an option for acute pain relief, an adjunct to physical therapy or rehabilitation, a post-surgical stabilization device, or simply a reminder to move carefully, the Aspen[®] LSO is effective, comfortable, and easy to use.

Sturdy Aluminum Rods

Reinforce spinal support and allow formation of patient-specific lordotic curves to minimize migration and optimize performance and fit.

Large Lumbar Window

Offers an unobstructed view of, and physical access to, the lower spine for custom curve assessment, better patient management, and true vented cooling.

SureSlot[™] Adjustment System

Allows for "gross" adjustment of circumferential size in a few seconds. When coupled with the LeverLatch[™] Tightener "fine" adjustment, the system ensures a precise fit for most patients by stocking only the two sizes.

Optional Cold Therapy

A reusable, soft gel cold pack slips into a thermal sleeve on a special pad to reduce pain or tenderness, and minimize swelling for better patient care, cooler wear, and comfort.



Performance You Can Count On

Aspen[®] LSO braces provide pain relief, stabilization, and comfort for performance that leads to better patient outcomes.